

endotine[®]
ribbon

ENGLISH

ENDOTINE Ribbon™ Instructions for Use

DESCRIPTION

The ENDOTINE Ribbon™ is a bioabsorbable fixation implant. This device is SINGLE USE ONLY and is supplied sterile.

INDICATIONS

The ENDOTINE Ribbon™ is intended for use in elevation and fixation of tissues in the temporal region, midface, lower face/jowl and the neck and to maintain fixation until healing has occurred. (Figure 1)

CONTRAINDICATIONS

- Situations where internal fixation is otherwise contraindicated (e.g. infection).
- Patients appearing to have very thin soft tissues of the face in which the implant may be visible or palpable.
- Any known allergy or foreign-body sensitivities to plastic biomaterial.

PACKAGING, LABELING AND STERILIZATION

The ENDOTINE Ribbon™ is pre-sterilized by gamma irradiation and intended for single use only. Do **NOT** resterilize this device or the insertion tool. The device and its insertion tool should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of puncture(s), tampering, water contamination or other damage.

STORAGE INSTRUCTIONS

Store at room temperature (15 to 24°C or 60 to 75°F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the label. Do not use if the temperature indicator dot has turned black; the indicator changes color if the device is subjected to temperatures that could damage it, even if the exposure is of short duration, as might occur during transportation.

INSTRUCTIONS FOR USE

- Figure 2 identifies the fixation region (2A), anchoring leash (2B), and integrated release tab (2C).
- The implant is packaged within the insertion tool (Figure 3), a protective cover (3A) through which the release tab (3B) protrudes.
- The dissection (temporal, midface, lower face/jowl, neck) is performed using the surgeon's typical approach (see appropriate Dissection Plane listed in Table A).

Indication	Dissection Plane	Time Orientation
temporal	Subcutaneous	Outward
midface	Subperiosteal	Outward
lower face/ jowl	Subcutaneous, times engaging SMAS fascia	Inward
neck	Superficial to platysmal muscle	Inward

Table A

- Insert the implant (tapered end first) through the incision into the desired position. If resistance is encountered, do not force the assembly into position, which may weaken the implant causing it to break. Remove the assembly and extend or enlarge your plane of dissection, with care to protect surrounding nerve structures, in order to provide adequate clearance for proper positioning of the device at the inferior recess of the dissected tissue.

- The release tab (Figure 3B) indicates the direction of the lines (see appropriate Time Orientation listed in Table A).
- When comfortable with the device position, retract the protective cover (Figure 4) by holding the release tab (4A) and removing the cover (4B).
- Apply digital pressure externally to engage the lines in the tissue (Figure 5). Suture may be placed in center holes to the distal fixation region for more secure fixation (optional, not necessary).
- To reposition the lines once they have been fully engaged in tissue, lift the skin over the lines with scissors, and disengage the lines with long pick-ups. Remove the device from the back of the cover (flat edge, not beveled), and redeploy.
- Place the device in the submandibular area within one to two finger breadths below the mandibular border.
- With the neck in the neutral position, engage and secure the device in the platysma by applying moderate (not extreme) tension to the anchoring leash exiting the incision.
- Once the desired position is reached, exercise caution with excessive device or tissue manipulation, endoscopic visualization, and while inserting or removing drains, in order to avoid puncture or damage the device. Avoid overlying implants to prevent palpability. Do not grasp the device with metal instruments, as this is a rubber shod Kelly clamp.
- After elevating the device to the desired position, place permanent, monofilament suture through the leash (Figure 6) and pass the suture through the soft tissue of adequate thickness. Secure the suture longitudinally in two locations, and vertically in two locations. Use caution if securing the Ribbon with permanent, braided suture. The braided suture may damage the device if it is tied too tightly.
- If a submental incision is made, the distal end of the Ribbon may be sutured for security.
- Trim the anchoring leash (Figure 6). It is desirable to bury the end of the leash either under adjacent, mobilized soft tissue or fascia to avoid possible erosion of the device through the skin. Close the incision using the surgeon's typical technique.
- If adjustment or removal of the implant is desired, one may remove the suture and extract the device under direct or endoscopic guidance.
- Postoperatively, the patient should be cautioned to avoid excessive neck rotation, rigorous physical activity such as aerobics or weight training and contact with the face and/or neck for the first two to three weeks post-procedure. The patient should be advised to wear a strong elastic garment for the first week post-procedure. After the first week, the patient should wear the elastic garment at night for an additional one to two weeks. Elevation of the head, neck, and shoulders while sleeping may reduce post-operative swelling.
- Removal of the device is required in the post-operative period, a small incision to release the device and remove the device using the most straightforward method.

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PRECAUTIONS

- Inform the patient of possible adverse effects including palpability of the device prior to surgical dissection, visibility of the device under the skin, and risk of infection inherent to the procedure (to include nerve injury, facial musculature dysfunction, infection, scarring, device erosion through the skin, asymmetry, loss of suspension/elevation, etc.). Some devices may require removal prior to their absorption due to discomfort, infection, reaction or other concerns. Some patients may desire temporary fillers to camouflage the borders of the implant during the healing phase.
- Transient local fluid accumulation and/or sinus formation may occur. Aspiration may yield traces of implant material and may result in healing of the sinuses without adverse effects.
- Use care around existing implanted devices.

WARNINGS

Do **NOT** RSTERILIZE the ENDOTINE Ribbon™ or its insertion tool. This device and cover are for SINGLE USE ONLY and MUST NEVER BE RE-USED. Resterilization may cause loss of function or failure of the device. Risks of re-use may include patient infection due to lack of sterility.

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