Annex to Motiva® Breast Implants Instructions for use.

ESTABLISHMENT LABS S.A.

Information for patients undergoing aesthetic breast augmentation
procedures

INFORMATION FOR THE PATIENT: BREAST IMPLANT SURGERY WITH MOTIVA IMPLANTS®

CAUTION

Only to be implanted in an appropriate medical environment by appropriately trained medical doctors who are qualified or accredited in accordance with national law.

Federal (USA) law restricts this device to be used by a certified plastic surgeon.

An extended version of this document is available through the website <https://www.motiva.health/ifu>.

1.INTENDED USE

MEDICAL PURPOSE (RECONSTRUCTION)

Motiva Implants® are intended to be used by surgeons under sterile conditions during breast augmentation (primary and revision).

2.INDICATIONS

NON-MEDICAL PURPOSE (AUGMENTATION)

Motiva Implants® are indicated in female patients of at least 18 years old for the following condition:

- To enhance the aesthetic appearance of patients dissatisfied with the shape/size of their breasts.
- To correct or improve the result of previous implant-based breast surgery.

3.INTENDED CONDITIONS FOR USE

Motiva Implants® must be used during breast implant procedures under sterile conditions, in compliance with good aseptic practices.

4.OVERVIEW

4.1. Alternative treatments are available for the conditions above, including external breast prostheses or padding or the transfer of other body tissues to enlarge breast size. Other synthetic filling materials (such as liquid silicone or other fillers) are not recommended and can provoke serious health problems.

4.2. The decision to have breast implants is a personal choice.

4.3. When choosing breast implant surgery, you should know that you may require additional procedures and further consultations with your surgeon. Silicone breast implants are not lifetime devices and are subject to wear and tear like any other implant device. Breast implantation might not be a one-time surgery. Your implant(s) may have to be removed or replaced, which may imply revision surgery. Many changes to your breasts following implantation are irreversible (cannot be undone). If you choose to have your implant(s) removed and not replaced, you may experience unsatisfactory aesthetic results, which can be permanent.

4.4. When you have your implants replaced, your risk of future complications increases compared to that associated with first-time breast implant surgery.

5.CONTRAINDICATIONS

The use of silicone breast implants is contraindicated in:

- Women with existing carcinoma of the breast, without mastectomy.
- Women with active infections.
- Women who are currently pregnant or nursing.
- Women with uncontrolled diabetes are clinically known to impact the wound-healing ability.
- Women who show tissue characteristics clinically incompatible with mammoplasties, such as tissue damage, compromised vascularity, or ulceration.
- Women with any condition – or treatment – determined by the surgeon to constitute an unjustifiable surgical risk (e.g., unstable cardiovascular disease, coagulopathies, chronic pulmonary problems, etc.).
- As a limitation, women with high levels of melanin have an increased risk of keloid and/or hypertrophic scarring.

6. RELEVANT TOPICS

6.1. Informed Consent

Establishment Labs relies on your surgeon to explain the existing risks and benefits of the implantation. The surgeon is also responsible for obtaining your formal informed consent to perform the surgical procedure.

To document a successful informed consent process, the surgeon must provide their “Informed Consent Document,” which must be signed by the surgeon, the patient, and a witness. This document will be part of the patient's medical file. Patients should be advised of all possible side effects and surgery-related complications when considering silicone gel-filled breast implants.

Summary of Safety and Clinical Performance (SSCP) - for the European Community, UK, and Switzerland:

The SSCP document could be consulted for additional and updated information on the Sterile Silicone Breast Implant Motiva Implants®, which will be available through the website <https://motiva.health/patients-motiva-experience/> and EUDAMED using the Basic UDI-DI.

7. APPROPRIATE MEDICAL ENVIRONMENT

This device should be used only in an aseptic environment in compliance with good surgical practices.

8. RISK AND BENEFIT ANALYSIS

8.1. EXPECTED CLINICAL BENEFITS

A benefit is expected when rebuilding or correcting the breast due to a clinical condition such as cancer, asymmetry, deformity, trauma or burns; since reconstruction offers many women the opportunity to restore the sense of balance and proportionality. Breast reconstruction eliminates the need for an external prosthesis or bra pads by creating a new breast that feels like a natural part of the patient's body and, it does not require harvesting tissue from another anatomical area. The procedure is customized to each patient's individual needs and situations. The information gathered from clinical studies confirms an improvement in patients' previous condition when used for breast reconstruction, as evidenced by the clinical outcome, and participant and surgeon's satisfaction appraisal.

No clinical benefits are expected for the aesthetic purpose of Motiva® Sterile Silicone Breast Implants. However, as a non-medical effect, patients dissatisfied with their breasts may expect an improvement in the aesthetic appearance.

8.2. RISKS OF BREAST SURGERY WITH SILICONE IMPLANTS

Motiva Implants® share the same intended use and mechanism of action to achieve the intended clinical functionality as other similar breast implants available in the market.

Motiva Implants® are not lifetime devices, and there is a possibility that patients will undergo implant removal(s), with or without replacement, throughout their lifetime. The longer you have your implants, the more likely it will be for you to have them removed or replaced, and the more likely you are to experience local complications and adverse outcomes. If you notice any abnormal changes in your breasts, you must see a doctor promptly. If you have silicone gel-filled breast implants, you will need to undergo periodic MR examinations to detect implant ruptures that do not cause symptoms (“silent ruptures”).

9. BREAST IMPLANT SURGERY SIDE EFFECTS

After breast implant surgery, patients might experience swelling, hardness, discomfort, itching, bruising, twinges, and pain over the first few weeks. The undesirable side effects identified are listed below.

9.1.Related to general anesthesia

General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and the facility where the surgery will occur to become aware of the tests, pre-surgical examinations, and length of time you need to be without food or your routine medications before the surgical procedure.

9.2.Residual risk

Residual risks	Description	Control measure
Inflammation/irritation	It is a foreign-body response demonstrated by redness, swelling, warmth, pain and/or loss of function.	Procedure-related adverse events associated with physiological factors; no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision. In addition, Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation process.
Gel fracture	It is defined as a fissure or crack in the implant's gel. It can occur due to exposing the implant to excessive compression forces applied to a small shell area during implant insertion, or due to the development of capsular contracture. It may result in device distortion.	Gel fracture can be detected by ultrasound or magnetic resonance imaging (MRI). Most gel fractures are undetectable by physical examination. Precautions regarding excessive force have been included.
Gel diffusion	Small quantities of silicone may diffuse/bleed through the shell of silicone gel filled implants.	BlueSeal® Technology provides a visual cue for quality control, enabling the confirmation of the barrier layer and its uniformity around all the shell areas.
Redness/erythema	Bleeding at the time of surgery can cause the skin to change color. It is likely temporary.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision.
Unsatisfactory results/ cosmetic defect	Unsatisfactory results such as stretch marks, visibility and dissatisfaction with the implant volume may occur.	Careful preoperative planning and surgical technique can minimize but not always prevent unsatisfactory results.
Pain	Most women undergoing breast implant surgery will experience post-operative breast and/or chest pain.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision.
Skin rash/silicone reaction/ allergy	Cutaneous risks with breast implants seem to be low. However, several reports have documented the presence of cutaneous hypersensitivity-like reactions to breast implants.	Topical and systemic medications may relieve the symptoms and lead to a successful resolution. In some cases, implant removal is required for complete symptom relief.
Swelling/edema	Breasts will feel tight after surgery and swelling will peak around three to four days later. However, you may experience some swelling for up to three months after surgery and it can take up to six months for the implants to fully settle in place.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision.
Galactorrhea	Is a rare event after implants-based breast surgery described in several cases in the literature, which consists of the spontaneous and usually bilateral milk production from the nipple unrelated to pregnancy, lactation or puerperium period.	Depending on the risk factors, symptoms, and imaging findings, treatment should be considered, which may be conservative (observation and surveillance), pharmacological and in some rare cases it must be treated with surgery.
Hematoma	It is a collection of blood within the space around the breast implant. Symptoms of hematomas generally include swelling, bruising and pain around the incision area. While most hematomas are small and will completely drain on their own and the blood re-absorbs into the body, those patients that are experiencing moderate to severe pain should undergo a follow-up visit.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision.
Bottoming-out	Inferior displacement of a breast implant that increases the distance between the nipple-areolar complex and the inframammary fold (IMF) after breast-implant surgery.	Appropriate surgical planning can mitigate the possible causes of bottoming-out. Recommendations include a careful and individual assessment of mammary tissues, careful implant selection, application of risk-minimizing surgical techniques, and adequate breast support after surgery.
Lactation difficulties	Women who undergo mastectomies and then have breast implant-reconstruction surgeries may not be able to breastfeed on the affected side due to loss of breast tissue and the glands that produce milk.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient; however, additional safety information is provided in the information for the patient for an informed decision. Precaution was included regarding periareolar incision that may considerably reduce the possibility of future breastfeeding.
Double capsule	It refers to the finding of two distinct capsular layers, separated by an inter-capsular space, around a breast implant. Clinical signs may vary from asymptomatic to firmness of the implant, discomfort, change in shape or position of the implant, and pain.	Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation and possibly, inflammation-related complications such as capsular contracture, double capsules, and seromas

Hypertrophic, keloid scarring and hyperpigmentation	Hypertrophic scars and keloids arise from some injury to the skin that goes deep enough to affect the dermal layer. Biologically, some people tend to be more susceptible to hypertrophic scars due to their genetic predisposition. Although uncommon, hyperpigmentation after surgery is a natural response to inflammation in some people.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Review contraindications during surgery planning regarding high melanin levels that have an increased risk of keloid and/or hypertrophic scarring.
Delayed wound healing	Some patients may experience a prolonged wound healing time. Smoking causes a decrease in the blood's oxygen levels, directly affecting the surgical wounds' healing process.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.
Breast mass/cyst	It is a fluid filled sac that develops within breast tissue. These sacs form when normal fluid producing glands in the breast enlarge or become blocked. Breast cysts are usually found on a self-breast exam. When they are small, they often go unnoticed or may instead be seen on a mammogram.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event.
Capsular contracture	Hypertrophic scar tissue investing in a foreign body or surgically implanted device, compromising the aesthetic outcome, resulting in pain, breast deformity, and often necessitating further operations.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Motiva® controlled surfaces have the advantage of less friction, with an improved performance in terms of capsular contracture and rupture rates, compared to conventional smooth and textured surfaces.
Infection	It can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. In addition, Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation process.
Changes in Nipple and Breast sensation	Breast surgery can result in an increased/decreased breast and/or nipple sensitivity. While some of these changes can be temporary, they can also be permanent and may affect the patient's sexual response or breastfeeding ability.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision.
Seroma	It is an accumulation of fluid that results from tissue inflammation. Often, seromas are reabsorbed by the body over several weeks, but needle drainage is sometimes needed to remove the fluid. Seroma symptoms most often appear a week to 10 days after surgery; the area may feel tender and swollen with a discrete lump and redness arising within a day or two.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation and possibly, inflammation-related complications such as capsular contracture, double capsules, and seromas
Symmastia	This is a relatively rare implant-displacement problem that occurs when the skin and muscle between the breasts over the sternum (breastbone) detaches and the two pockets of tissue that hold the breast implants come together to form one pocket. This allows the implants to come together in the middle, creating the appearance of a "uniboob" and sometimes causing discomfort or pain. It's often difficult to correct symmastia and it may require more than one surgery.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision.
Breast Atrophy/Chest Wall Deformity	The implant's pressure may cause the breast tissue to thin and shrink (with increased device visibility and palpability), potentially leading to chest wall deformity	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision.
Rupture	Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation but is more likely to occur the longer the implant is in place. The integrity of breast implants (and detection of gel fractures and/or silent ruptures) can be evaluated through several techniques. High-resolution ultrasound (HRUS) is widely accepted by healthcare providers and patients for rupture diagnosis. Additionally, the USFDA recommends magnetic resonance imaging (MRI) surveillance with the first MRI performed three years postoperatively and subsequent MRIs performed every two years after that.	Motiva® controlled surfaces have the advantage of less friction, with an improved performance in terms of capsular contracture and rupture rates, compared to conventional smooth and textured surfaces.

Extrusion	Breast-implant extrusion, or breast implant exposure, occurs when the breast skin and tissues that are holding the implant fail causing the implant to protrude through the skin and become exposed. Breast-implant extrusion calls for surgery and removal of the implant	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.
Surgical Wound Dehiscence (SWD)	In a superficial wound, the risks of Surgical Wound Dehiscence (SWD) increase, causing a separation of the margins of a closed surgical incision that has been made in the skin with or without exposure or protrusion of underlying tissue, organs, or implants. The separation may occur in one or multiple regions, involve the entire length of the incision, and affect one or more layers of tissue.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.
Rotation	Anterior/posterior rotation, also called flipping, has been observed more frequently with cohesive gel implants. The flat base of the implant is positioned anteriorly, deforming the breast of the patient. In some cases, revision surgery may be necessary to reduce pocket dimensions. Diamond projection of Ergonomix2® has a symmetrical shape with the same apex on both sides, a designed geometry that could reduce visible profile divergences in the event of flipping.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Surgical technique can minimize the possibility of this harm.
Calcification/calcium deposits	Calcium deposits can form in the scar tissue surrounding the breast implants and cause pain and firmness. Calcium deposits are visible in mammography. As these deposits must be differentiated from calcium deposits that are a sign of breast cancer, additional surgery may be necessary to remove and examine calcifications.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision.
Malposition	It is defined as an incorrect placement during surgery or its shifting from its original position. Malposition has been a frequent event reported due to its multifactorial causes and it can be expected during the lifetime of the device. The implant's shifting can be produced by trauma, capsular contracture, gravity, or initially wrong placement with the device. The risk associated with this event is dissatisfaction with aesthetic outcomes.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Surgical technique can minimize the possibility of this harm.
Interference with mammography	You should be instructed to inform your examiners about the presence, type and placement of the implants and request diagnostic mammography rather than screening mammography. Breast implants may complicate the interpretation of mammographic images by obscuring underlying breast tissue and/or compressing overlying tissue.	Communication with radiologist or technician regarding the presence of implants to use appropriate techniques. Information related to this known adverse event has been included in the patient information and Directions for Use.
Asymmetry	Postoperative aesthetic difference in the two breasts produced by different factors, such as a fall of the fold, a high implant, or a rotation of the implant. It is possible that after breast implant surgery, small deformities in the wall of the thorax or a morphologic mammary disorder may become much more evident.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on morphological characteristics.
Necrosis	It refers to the formation of dead tissue around the implant. This may prevent wound healing and require surgical correction and/or implant removal. Factors associated with necrosis include infection, steroids in the surgical pocket, chemotherapy/radiation, and excessive heat or cold therapy	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.
Rippling/wrinkling	It is the cutaneous manifestation, visible or palpable, of the implant ripples and edge that are typically most apparent when the patient bends forward. In situations where the implant's soft tissue coverage is insufficient, these harmful effects become more apparent.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Surgical technique can minimize the possibility of this harm. Precaution was included regarding wrinkles, it is recommended to run a finger around the implant before closing to ensure the implant is flat.
Ptoxis/drooping/waterfall effect	The "waterfall effect" is a descriptive term to indicate sliding ptosis of parenchymal breast tissue over a fixed or encapsulated implant. Certain breast implants are more prone to contribute to this problem as are implants placed in submuscular pockets that ride high, especially in women with anatomical musculoskeletal variance or asymmetry.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.

EN		
Reoperation (Explanation)	Rupture, unacceptable cosmetic outcomes (i.e., dimpling, wrinkling, and other potentially permanent cosmetic changes of the breast), and other complications may require additional surgeries to the patient's breasts. Implants are not lifetime devices, and there is a possibility that patients will undergo implant removal(s), with or without replacement, throughout their life. When implants are explanted without replacement, changes to the patient's breasts may be irreversible.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision.
Breast cancer	The Food and Drug Administration (FDA) identified a possible link between breast implants and anaplastic large cell lymphoma (ALCL), an uncommon cancer of the immune system. When this cancer occurs, it is called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL).	Review contraindications during surgery planning. Information related to this known adverse event has been included in the patient information and Directions for Use.
Visibility/palpability	Visibility and palpability of the implant can occur because the envelope is thin due to several causes, including excessive volume. The implant's content not being cohesive, or previous surgery and cutaneous aging.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event based on the quality of the tissue.
Lymphadenopathy	Silicone-induced lymphadenopathy is a well-known, rare complication of implant insertion. It is a disease of the lymph nodes (i.e., the small, round structures that operate as part of the body's immune system). They become abnormal in size or consistency (most commonly producing swollen or enlarged lymph nodes). After implant insertion, axillary lymphadenopathy causes are multifactorial with possible factors, including granulomatous reaction, inflammation and/or malignancy. Breast-implant rupture and/or leakage through an intake surface can cause fibrosis and granulomatous reactions, which in turn can result in a contracture or regional lymphadenopathy that sometimes mimics malignancy.	Device and procedure related adverse event associated with physiological factors or foreign material, the information has been included in the patient information.
Connective tissue disease (CTD)	The evidence has not detected any affirming association between silicone gel-filled breast implants and connective tissue disease or any other autoimmune diseases. However, a spectrum of autoimmune symptoms such as fatigue, myalgias/arthralgias, dry eyes/mouth, rash, "brain fog," neurological disturbances (e.g., visual symptoms or alterations in sensation, muscle strength, walking, balance, thinking or memory) or diseases (e.g., multiple sclerosis), have been reported by some patients with both saline and silicone gel-filled breast implants in the last review of the medical report databases from FDA.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event.
Interference with Magnetic Resonance Imaging (MRI): Missing an anatomical area of interest due to artifact leading to misdiagnosis or missed diagnosis / delayed treatment (non-ferrite microtransponder-Zen)	Motiva Implants® with a microtransponder are classified MRI-conditional. The microtransponder, during an MRI scan, can create an image immediately around the microtransponder known as a susceptibility artifact. This interference may cause visualization challenges when reading tissue structures and can present as shadows, distortions, or voids within the imaging field. When a dual-modality approach (MRI and US) is offered to the patients who belong to the breast cancer high-risk group, on average one out of 17,892 exams would likely miss a patient with cancer recurrence. As such, the radiologist and/or radiology technologist should be informed of the presence of the microtransponder, notifying them that it is embedded in the posterior part of the device near the patch of the breast implant. This will allow the professional to anticipate the presence and accurately manage an imaging artifact with certain MRI techniques.	MRI Safety Indications are available for surgeons and radiologists. Other imaging modalities that may be used in addition to MRI and are recommended to complement visualization and improve overall diagnosis, such as: Breast ultrasound, Digital tomosynthesis, Automated ultrasound - Scintimammography, Dual-energy contrast subtraction digital mammography, Dual-energy enhanced mammography.
Breast implant illness (BII)	Patients and clinicians use the term "breast implant illness" or "BII" to describe a variety of systemic symptoms that have been reported by some women following reconstruction or augmentation with breast implants. The possible symptoms related include fatigue, arthralgia, muscle pain, memory loss, difficulty concentrating, rash, dry eyes, brain fog and visual disturbance, and autoimmune disorders. BII is not recognized as a formal medical diagnosis, and there are no specific tests or recognized criteria to define or characterize it. However, studies have shown that silicone is a biologically inert material and unlikely to be the cause of these symptoms.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. The review of contraindication during preoperative planning could minimize this adverse event.
Breast implant associated squamous cell carcinoma (BIA-SCC)	There have been literature reports of SCC in the capsule formed around breast implants for both textured and smooth breast implants and for both saline and silicone breast implants. The cases were reported in patients who were diagnosed 15–42 years after having breast implants placed and have all presented with abnormal signs and symptoms, including swelling, lumps, or skin changes. The incidence of BIA-SCC is very rare; and according to the case reports, it behaves similarly to BIA-ALCL, with local capsular growth followed by systemic metastasis.	Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation and possibly, inflammation-related complications such as capsular contracture, double capsules, and seromas

EN	
Breast implant associated anaplastic large cell lymphoma (BIA-ALCL)	It is not breast cancer - it is a type of non-Hodgkin's lymphoma (cancer of the immune system). In most cases, BIA-ALCL is found in the scar tissue and fluid near the implant, but in some cases, it can spread throughout the body. An individual's risk of developing BIA-ALCL is considered to be low; however, this cancer is serious and can lead to death, especially if not treated promptly. The main symptoms of BIA-ALCL are persistent swelling, a mass, or pain in the breast implant area. Most cases of BIA-ALCL occur years after breast implant placement and present with symptoms or changes around the breast implant. It has been reported that BIA-ALCL occurs more frequently following implantation of breast implants with textured surfaces rather than those with smooth surfaces. In cases of ALCL, evidence of fluid collection around the breast implant (seroma) is often observed during evaluation by a health care provider. Some patient reports indicated that a lump under the skin or capsular contracture (thick and noticeable scar capsule around the implant) were present. Procedure-related adverse event associated with physiological factors, no specific control measures can be taken by the patient, however, additional safety information is provided in the information for the patient for an informed decision. Motiva® controlled surfaces were designed to improve compatibility between implants and tissues, minimizing inflammation and possibly, inflammation-related complications such as capsular contracture, double capsules, and seromas In most patients, it is treated successfully with surgery to remove the implant and surrounding scar tissue, and in some patients, also treatment with chemotherapy and radiation therapy. Most cases of BIA-ALCL are treated by removal of the implant and the capsule surrounding the implant, and some cases have been treated by chemotherapy and radiation.
	10. WHAT TO EXPECT FOLLOWING SURGERY (GENERAL POST-OPERATIVE CARE) During your recovery process, you may experience the following: - You will likely feel tired and sore for several days following the operation. - Your breasts may remain swollen and sensitive to physical contact during expansion. - You may experience a feeling of tightness in the breast area as the skin adjusts to the new breast size, which may increase after each expansion. - After surgery, some patients report a higher body temperature. You should notify your surgeon to determine whether this is a normal response or if you may need to seek medical attention. It is important to notify your surgeon if you have a fever and any of the following symptoms: - Nausea - Vomiting - Unexplained worsening of pain, disorientation, or confusion - Purulent or serous drainage from your incision - Redness and/or bruising near your incision. - Shortness of breath - After surgery, you will sometimes be left with one or two drains, which prevent the accumulation of bloody or serous fluid in the tissues. Both may remain in place for 5-15 days or until your surgeon decides to remove them. Take all post-operative antibiotics and other medications as prescribed by your surgeon. - Sleep or rest with your head slightly elevated, avoiding lateral positions. - Keep your arms close to your body and avoid lifting weights until your surgeon allows it. - You should not drive for at least two weeks or perform physical exercise or gymnastics for at least four weeks or until authorized by your surgeon. - Do not use deodorant or other cosmetic products until you have worn your bra for some time. This will prevent deodorant from entering your incision and possibly causing an infection. Do not use deodorant on the side of your surgery if there is any break in the skin. - Do not expose your breasts directly to sunlight until your surgeon clears you. Sun exposure may cause irreversible skin discoloration or increase inflammation of the reconstructed breast. - If you are instructed to use an elastic bandage or surgical bra, only remove it for bathing. It should be used every day for 24 hours for at least six weeks or as long as your surgeon instructs you. - The bra or bandage should NOT be tight or unworned to allow proper breathing and not cause pain due to pressure, nor should it be too loose or baggy. - Pregnancy and nursing after breast implant surgery may cause breast tissue and muscle changes that could lead to ptosis (drooping) and flapping. - If any issues are present, you should report them immediately and possibly have an MRI evaluation to screen for rupture. 11. STERILE PRODUCT Motiva Implants® are sterilized during manufacturing by using a dry heat sterilization method. An implant is intended to be used only in one patient for a single procedure and supplied in a sealed, double sterile barrier primary package. 12. FOLLOW-UP EXAMINATIONS 12.1. Symptomatic Rupture Symptoms associated with rupture may include hard knots or lumps surrounding the implant, loss of size, pain, tingling, swelling, numbness, burning, or hardening of the breast area. If you notice any of these changes, consult your plastic surgeon to see that she/he can examine your implant(s) for rupture and determine whether you need to have an MR examination to find out if your symptoms are due to implant rupture. If rupture has occurred, you should have your implant removed and/or replaced. 13. IMPORTANT INFORMATION FOR PATIENTS UNDERGOING MAGNETIC RESONANCE IMAGING (MRI) In some cases, performing an MRI may be necessary. It is up to your surgeon and/or treating physician to determine and complete the MRI testing as required by your clinical condition. Patients with Motiva Implants® may safely undergo MRI testing under certain conditions. If you have been instructed to have an MRI, provide your MRI technician with the important safety information described in Appendix 1 of this document. Your MRI technician should review this important safety information to determine if it is safe for you to undergo MRI testing while your breast implant is in place. You will be monitored continuously throughout the MRI procedure using visual and audio means (e.g., intercom system). You will be instructed to alert the MRI technician of any unusual sensations or problems so that, if necessary, he or she can immediately terminate the procedure. - Heat Sensation- Sterile Silicone Breast Implants Motiva Implants® with microtransponder, during an MRI scan under defined conditions, can produce a minimal heat sensation 14. ADDITIONAL INFORMATION 14.1. Life Expectancy Motiva Implants® are not lifetime devices. Based on a vast literature review from clinical data obtained from equivalent devices, a survival rate in ten years of 80% has been established as acceptance criteria for the Motiva Implants® lifespan. 14.2. Device Traceability Motiva Implants® are subject to device tracking via the MotivaImagene® registration system. You can register your implants at https://motiva.health/motivaimagene/. If you have difficulty registering your implants, you can contact Establishment Labs to receive assistance. 14.3. Product Evaluation Establishment Labs requires that any complications resulting from the use of Motiva Implants® be reported immediately to your doctor. Your doctor must fill out all the necessary information using the Motiva Implants® Complaint Form available at the following webpage: www.motiva.health/support.

14.4. Patient Implant Card
You must have a record of your surgical procedure in case of future consultations or additional surgeries. Each implant comes with a Patient Implant Card, which must be given to you by your surgeon for personal reference. Besides the information stated on the Patient Record Label (which should come affixed to the back of the card), the Patient Implant Card also includes your name (patient identification), the position of the implant, date of implantation (surgery date), and the name of the treating surgeon (healthcare center or doctor). This implant card is for patients' permanent records and should always be kept safely.

NON-HARMONIZED SYMBOLS IN PRODUCT LABELING

Below, you can find symbols for your reference.

	Patient Identification		Imaging Studies Conditional		Implant volume
	Healthcare center or doctor		Manufacturer		Unique Device Identifier
	Date		MR conditional, the device can be imaged safely under the tested specifications described in Directions for Use.		Medical device
	Position of the implant		Catalog number		Serial Number

MAIN STANDARDS AND COMMON SPECIFICATIONS APPLIED

Code	Description
(BS EN) ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes.
(BS EN) ISO 14971	Medical devices - Application of risk management to medical devices.
(BS EN) ISO 14607	Non-active surgical implants. Mammary implants. Particular requirements
(BS EN) ISO 10993	Biological evaluation of medical devices - Parts 1, 3, 4, 5, 6, 10, 11 and 18.
(BS EN) ISO 14630	Non-active surgical implants. General requirements.
(BS EN) ISO 11607-1	Packaging for terminally sterilized medical devices — Requirements for materials, sterile barrier systems and packaging systems
(BS EN) ISO 11607-2	Packaging for terminally sterilized medical devices — Validation requirements for forming, sealing and assembly processes
(BS EN) ISO 20857	Sterilization of health care products — Dry heat — Requirements for the development, validation and routine control of a sterilization process for medical devices
(BS EN) ISO 11737	Sterilization of health care products - Microbiological methods- Parts 1 and 2.
(BS EN) ISO 62366-1	Medical devices-Application of usability engineering to medical devices
(BS EN) ISO 15223-1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied
(BS EN) ISO 20417	Medical devices — Information to be supplied by the manufacturer
Regulation(EU) 2022/2346	Common specifications for the groups of products without an intended medical purpose listed in Annex XVI to Regulation (EU) 2017/745 of the European Parliament
Note: Standards and common specifications are applied in the last version released or published	

15. REPORTING AND ADDITIONAL INFORMATION

If you need additional information related to Motiva Implants®, do not hesitate to contact us. If any serious incident occurs, go immediately to your surgeon and report the event to the closest Establishment Labs office:

LEGAL MANUFACTURER

Establishment Labs S.A.
Coyol Free Zone & Business Park Building
4th Street, Building B-15, Alajuela, Costa Rica.
Zip code: 20113
Phone: +506 2434-2400
www.motiva.health/support
www.motiva.health

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Applicable to patients in EU Members States, UK, and Switzerland:

Any serious incident related to Motiva Implants® devices should be reported to Establishment Labs and the corresponding competent authority in which the patient is established, if applicable according to local regulations.

APPENDIX 1. MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

MRI Safety Information

A person implanted with the Motiva Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury. Additionally, MRI safety information can be readily accessible in the Motiva website.

Device Name	Motiva Implants® with microtransponder
Static Magnetic Field Strength (B ₀)	1.5-Tesla (64-MHz) and 3-Tesla (128-MHz) only
Maximum Spatial Field Gradient	40- T/m (4,000-gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	There are no Transmit Coil restrictions
Operating Mode	Normal Operating Mode

Maximum Whole-Body Specific Absorption Rate (SAR)	2-W/kg in the Normal Operating Mode
Maximum Head SAR	3.2 W/kg in the Normal Operating Mode
Scan Duration	2-W/kg whole-body average SAR for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode 3.2-W/kg head average SAR for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode
MR Image Artifact	The presence of this implant may produce an image artifact.
Additional Precautions	Under the scan defined conditions, the Motiva Implants® are expected to produce a maximum temperature rise of 1.5°C after 15-minutes of continuous scanning (i.e., per pulse sequence). The Motiva Implants® will not present an additional risk or hazard to a patient in the 3-Tesla MR environment concerning translational attraction or migration. Do not perform the MR if the patient is sedated, anesthetized, confused, or otherwise unable to communicate with the MR system operator. If information about a specific parameter is not included, please contact the manufacturer.

APPENDIX 2. DETAILS OF THE VOLUME AND SIZE OF APPLICABLE DEVICES.

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Copyright 2025 Establishment Labs S.A.