

Liposuction Cannulas

MICROAIRE®

ASPIRATION

Mercedes

- 360° of extraction
- Moderately aggressive
- The most popular cannula style

Double Mercedes

- 360° of extraction
- Aggressive
- Good for fibrous tissue and large-volume extraction

Flared Mercedes

- 360° of extraction
- Aggressive
- Ideal for debulking, dissection, and pre-tunneling
- Good for breaking up fat without suction, prior to sculpting and extraction

Bent Flared Mercedes/Mercedes

- Effective for pre-tunneling
- Good for debulking and extraction

Bent Flared Mercedes

Bent Mercedes

Tri-Port II

- 270° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Mirrored Tri-Port II

- 270° of extraction
- Aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Tri-Port III

- 90° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Helixed Tri-Port III

- 270° of extraction
- Aggressive
- Designed for high-fibrous areas
- Ideal for gynecomastia cases

Single-Use	Multi-Use	Diameter	Length	Port Length
PAL-504LL	PAL-R504LL	5mm	30cm	10mm
PAL-504LS	PAL-R504LS	5mm	22cm	10mm
PAL-404LL	PAL-R404LL	4mm	30cm	10mm
PAL-404LS	PAL-R404LS	4mm	22cm	10mm
PAL-304LL	PAL-R304LL	3mm	30cm	8mm
PAL-304LS	PAL-R304LS	3mm	22cm	8mm
PAL-304LM	PAL-R304LM	3mm	15cm	8mm
PAL-244LM	PAL-R244LM	2.4mm	15cm	7.5mm
—	PAL-R404XL	4mm	40cm	10mm
PAL-505LL	PAL-R505LL	5mm	30cm	10mm
PAL-505LS	PAL-R505LS	5mm	22cm	10mm
PAL-405LL	PAL-R405LL	4mm	30cm	10mm
PAL-405LS	PAL-R405LS	4mm	22cm	10mm
PAL-508LL	PAL-R508LL	5mm	30cm	10mm
PAL-508LS	PAL-R508LS	5mm	22cm	10mm
PAL-408LL	PAL-R408LL	4mm	30cm	10mm
PAL-408LS	PAL-R408LS	4mm	22cm	10mm
PAL-308LL	PAL-R308LL	3mm	30cm	8mm
PAL-308LS	PAL-R308LS	3mm	22cm	8mm
PAL-508LLB	PAL-R508LLB	5mm	30cm	10mm
PAL-408LLB	PAL-R408LLB	4mm	30cm	10mm
PAL-404LLB	PAL-R404LLB	4mm	30cm	10mm
PAL-304LLB	PAL-R304LLB	3mm	30cm	8mm
PAL-402LL	PAL-R402LL	4mm	30cm	10mm
PAL-402LS	PAL-R402LS	4mm	22cm	10mm
PAL-302LL	PAL-R302LL	3mm	30cm	8mm
PAL-302LS	PAL-R302LS	3mm	22cm	8mm
PAL-506LL	PAL-R506LL	5mm	30cm	10mm
PAL-506LS	PAL-R506LS	5mm	22cm	10mm
PAL-406LL	PAL-R406LL	4mm	30cm	10mm
PAL-406LS	PAL-R406LS	4mm	22cm	10mm
PAL-403LL	PAL-R403LL	4mm	30cm	10mm
PAL-403LS	PAL-R403LS	4mm	22cm	10mm
PAL-303LL	PAL-R303LL	3mm	30cm	8mm
PAL-303LS	PAL-R303LS	3mm	22cm	8mm
PAL-507LL	—	5mm	30cm	10mm
PAL-507LS	—	5mm	22cm	10mm
PAL-407LL	PAL-R407LL	4mm	30cm	10mm
PAL-407LS	PAL-R407LS	4mm	22cm	10mm

				Single-Use	Multi-Use	Diameter	Length	Port Length
Del Vecchio Track-12  • 270° of extraction • Aggressive • Designed for harvesting small lobules for grafting				PAL-4012LL	—	4mm	30cm	6mm
				PAL-3012LS	—	3mm	22cm	6mm
				—	PAL-R4012XL	4mm	40cm	6mm
Multi-Hole  • 360° of extraction • Non-aggressive • Produces small lobules for small-volume grafting • Also ideal for fine sculpting and feathering				PAL-30MLL	PAL-R30MLL	3mm	30cm	1.9mm
				PAL-30MLS	PAL-R30MLS	3mm	22cm	1.9mm
				PAL-30MLM	PAL-R30MLM	3mm	15cm	1.9mm
Spatula  • 90° of extraction • Short, non-aggressive cannula with a blunt tip • Ideal for delicate areas such as the neck and back	Three-Hole	---		PAL-30S3	—	3mm	15cm	1.5mm
	Two-Hole	---		PAL-30S2	—	3mm	15cm	1.5mm
	Single-Hole			PAL-24S1	—	2.4mm	15cm	3mm
Single-Port  • 90° of extraction • Non-aggressive • Good for sensitive areas such as the neck and back • Ports can be directed away from the skin to prevent surface irregularities				PAL-300LL	—	3mm	30cm	8mm
				PAL-300LS	—	3mm	22cm	8mm
				PAL-300SS	—	3mm	22cm	6mm
				PAL-240LS	—	2.4mm	22cm	7.5mm
				PAL-240SS	—	2.4mm	22cm	5mm
				PAL-240LM	—	2.4mm	15cm	7.5mm
Turbo  • Allows air to vent through the tube during liposuction for better flow rates • Not recommended for harvesting fat to be re-injected	Double Mercedes	---		PAL-505LLT	—	5mm	30cm	10mm
				PAL-405LLT	—	4mm	30cm	10mm
	Helixed Tri-Port III	---		PAL-507LLT	—	5mm	30cm	10mm
	Mercedes			PAL-404LLT	—	4mm	30cm	10mm
				PAL-404LST	—	4mm	22cm	10mm
HD  • 270° of extraction • Moderately aggressive • Ideal for sculpting of superficial fat layer	Left Bent	---		PAL-40HLLB	—	4mm	30cm	10mm
	Straight	---		PAL-40HLL	—	4mm	30cm	10mm
	Right Bent			PAL-40HLRB	—	4mm	30cm	10mm

INFILTRATION

Multi-Hole  • 360° of infiltration				PAL-INF-30MLL	—	3mm	30cm	1.9mm
				PAL-INF-30MLS	—	3mm	22cm	1.9mm