

PAL[®] Infiltration System



PAL[®]
BY MICROAIRE

PAL® Infiltration System

The PAL® Infiltration Pump integrates with the PAL® System to activate the peristaltic pump from the PAL Handpiece. Custom electronic controls ensure consistent flow at low speeds and a cooling fan prevents overheating of the motor from heavy use. PAL® Infiltration Tubing eliminates obstructions with a large-bore design and smooth connections.



PAL® Infiltration Pump
INF-PUMP



Pump Type	Peristaltic
Average Max Flow Rate	1,155 mL/min
Weight	9.3 lbs (4.2 kg)
Dimensions	11.2" x 10.1" x 5.1" (29 x 26 x 13 cm)

System Accessories




Infiltration Pump to 5020
Connector Cable
INF-CBL-5020



Infiltration Foot Switch
INF-FOOT-1



PAL® Infiltration Tubing
PAL-INF-1600



PAL® Single-Use Infiltration Cannulas
PAL-INF-30MLS
PAL-INF-30MLL

A versatile system that expands
the capabilities of PAL®



Ordering Information

PAL® Infiltration System		Qty.
INF-PUMP	PAL® Infiltration Pump	EA
INF-CBL-5020	Infiltration Pump to 5020 Connector Cable	EA
INF-FOOT-1	Infiltration Foot Switch	EA
PAL-INF-1600	PAL® Infiltration Tubing	5PK
PAL-INF-30MLS	Single-Use Multi-Hole Infiltration Cannula, 3mm diameter, 22cm length	EA
PAL-INF-30MLL	Single-Use Multi-Hole Infiltration Cannula, 3mm diameter, 30cm length	EA

To order, call toll free **800·722·0822**

Learn more about our products, including the PAL® Infiltration System, at www.microaire.com.

MICROAIRE®

MicroAire Surgical Instruments
3590 Grand Forks Boulevard
Charlottesville, VA 22911 USA
800 722 0822 or 434 975 8000
www.microaire.com

©2025 MicroAire Surgical Instruments, LLC
All rights reserved. LIT-INF-PUMP Rev C



A Marmon |
Berkshire Hathaway
Company